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Design of an Adaptive Optical System

原文：Flexible Optical B.V.,The Netherlands（OKO Tech）

翻译：灵动智能光学，152 5716 7514 , info@inspiroptics.com

1. Phase Conjugation（相位共轭）

The optical scheme of a phase-conjugated AO system should satisfy the following conditions.

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（1）.The optics should re-image the plane of the mirror to the plane of the Hartmann mask (or microlens array) with proper scaling. An example of such an optical conjugation scheme is shown in Figure 1.

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（2）. In a similar way, optical conjugation should be provided between the DM and the source of aberrations.

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（3）. In a similar way, optical conjugation should be provided between the system pupil, and the tip-tilt corrector, if any.

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In general, phase-conjugated AO contains three planes conjugated to the pupil: the DM, the WF sensor and the tip-tilt corrector. Each conjugation requires a telescope, and each telescope introduces additional static aberration and scattering.

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Typical measurement scheme for testing of reflective optics with Shack-Hartmann wavefront sensor. Here CL is a collimator; L1 and L2 are lenses with focal distances f1 and f2, respectively; WFS is the wavefront sensor.

To reduce the number of conjugated planes, OKO fabricates MMDM15-17-tip-Tilt—special deformable mirror mounted on a tip-tilt stage.

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An important remark should be done regarding the use of micromachined membrane deformable mirrors (MMDM). As the membrane can only be pulled by electrostatic forces but cannot be pushed, it should be biased to provide deformation in both directions. It means that one will get much better results when optimizing the mirror with respect to a sphere, than when optimizing it with respect to a flat surface. This sphere should be taken into account when the mirror is incorporated into the optical setup. The optimum radius of the bias curvature can be evaluated assuming the mirror flexure to be in the middle of the range.

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Scheme of typical adaptive optics setup.

The typical setup for functional feedback loop is shown in the Figure 2. A high-quality flat or spherical mirror can be used as a reference. A concave spherical mirror is preferable for the use with MMDM due to its biased operation, and a flat mirror for PDM systems.

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2. FrontSurfer in the AO mode (用 FrontSurfer 软件进行 AO 闭环校正)

Closed-loop correction mode is implemented in the version of“FrontSurfer”for deformable mirrors.

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Before starting closed-loop correction with a deformable mirror (DM), the mirror should be properly configured. “FrontSurfer” supports standard OKO’s ISA, PCI, USB and ethernet interface modules for deformable mirrors. Configuration files can be supplied for standard OKO’s DM types. It is possible to integrate “FrontSurfer” with a non-OKO DM interface by developing a custom plugin according to our specification.

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2.1 Calibration of the AO (AO 系统校准)

After the mirror is coupled with the wavefront sensor, it should be calibrated.

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The calibration consists in the measurement of all Shack-Hartmann sensor responses (displacements of spots) corresponding to the actions of the actuators. These responses form the influence matrix of the mirror; we shall denote it as A. Calibration should be performed after any change in the optical setup.

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Before calibration, one needs to make sure that:
1. The DM is optically conjugated to the WF sensor;
2. The pupil of the system is optically conjugated to the DM.

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The first condition is extremely important for the stability of the feedback. Bad alignment of conjugated planes causes geometrical distortions and movement of the beam on the WF sensor, resulting in an unstable behavior of the AO system. The second condition is important, as it defines the final quality of the AO correction. Conjugation between the DM and the WF sensor can be checked by introducing tilts in the mirror plane. These tilts should not cause any movement of the beam in the sensor input. In the HS response one should observe only movement of spots. If the conjugation is correct, no new spots, caused by the beam movement on the sensor input, should appear with the DM tilts.

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After the pupil images are aligned and conjugated, the reference spot pattern, corresponding to the ideal performance of the AO system, should be obtained. The method of obtaining such an image depends on the AO system structure. Usually it is registered by injecting a diffraction limited reference beam into the the AO system in such a way, that all internal aberrations of the system are accounted for in the reference image. In the Fig. 2 it is achieved by blocking the DM arm and registering the spot pattern obtained by reflection on the reference mirror. In practice,obtaining the reference image that takes into account all the internal aberrations of the AO

system can happen to be one of the most important and cumbersome operations in the whole process of the AO system alignment and calibration.

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After calibration, the influence matrix **A** is subjected to singular value decomposition (SVD) $\mathbf{A} = \mathbf{USV}^{\mathbf{T}}$ [1]. The columns of the matrix **U** make up an orthonormal set of the mirror deformations (modes), and the values of the diagonal matrix **S** represent the gains of these modes. Discarding modes with small singular values improves controllability of the system [2], however can lead to the accumulation of errors that are correctable by the mirror, but excluded from the set of SVD modes.

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The “FrontSurfer” software allows to visualize the modes and singular values, and to limit the number of modes employed for the WF correction.

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Singular values of a 37-channel MMDM.

Experimental singular values for a 37-channel micromachined membrane deformable mirror are given in Figure 3; first 20 SVD modes are shown in Figure 4.

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First 20 SVD modes of a 37-channel MMDM.

“FrontSurfer” perform wavefront correction in a series of iterations. If the residual aberration ϕ_n at the nth iteration corresponds to the set of actuator signals $\mathbf{X_n}$ then the actuator signals at the next step $\mathbf{X_{n+1}}$ will be given by

$$\mathbf{X_{n+1} = X_n - gA^{-1}}\phi_n,$$

where g is the feedback coefficient with value in the range (0..1], $\mathbf{A^{-1}}$ is the pseudo-inverse of A given by $\mathbf{A^{-1} = VS^{-1}U^T}$.

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The target wavefront is defined as a combination of Zernike polynomials, and can be manually controlled during the correction using arrow keys of the computer keyboard.’ “Left” and “Right” switch between Zernike modes while “Up” and “Down” control the amplitude of the chosen Zernike mode. Examples of Zernike aberrations with the maximum possible amplitude for which rms error was less than $\lambda/14$ (Maréchal criterion) generated in this way by 50-mm PDM with 39 actuators are shown in Figures 5-20.

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Initial aberration of the system: reconstructed wavefront, simulated interferogram, simulated far field, and registered interferogram.

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Flatten mirror: reconstructed wavefront, simulated interferogram,simulated far field, and registered interferogram.

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Tilt (Zernike term Z[1,-1]) of amplitude 7.6μm generated : reconstructed wavefront,simulated interferogram, simulated far field, and registered interferogram.

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Defocus (Zernike term Z[2,0]) of amplitude 2.2μm generated: reconstructed wavefront,simulated interferogram, simulated far field, and registered interferogram.

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Astigmatism (Zernike term Z[2,2]) of amplitude 4.8μm generated: reconstructed wavefront, simulated interferogram, simulated far field, and registered interferogram.

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Astigmatism (Zernike term Z[2,-2]) of amplitude −5.4μm generated: reconstructed wavefront, simulated interferogram, simulated far field, and registered interferogram.

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Coma (Zernike term Z[3,1]) of amplitude 1.2μm generated: reconstructed wavefront,simulated interferogram, simulated far field, and registered interferogram.

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Coma (Zernike term Z[3,-1]) of amplitude 1.3μm generated: reconstructed wavefront,simulated interferogram, simulated far field, and registered interferogram.

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Trefoil (Zernike term Z[3,3]) of amplitude 3.9μm generated: reconstructed wavefront, simulated interferogram, simulated far field, and registered interferogram.

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Trefoil (Zernike term Z[3,-3]) of amplitude 3.9μm generated: reconstructed wavefront, simulated interferogram, simulated far field, and registered interferogram.

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Spherical aberration (Zernike term Z[4,0]) of amplitude 0.7μm generated: reconstructed wavefront, simulated interferogram, simulated far field, and registered interferogram.

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Zernike term Z[4,2] of amplitude 0.6μm generated: reconstructed wavefront, simulated interferogram, simulated far field, and registered interferogram.

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Zernike term Z[4,-2] of amplitude 0.7μm generated: reconstructed wavefront, simulated interferogram, simulated far field, and registered interferogram.

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Zernike term Z[4,4] of amplitude 1.9μm generated: reconstructed wavefront, simulated interferogram, simulated far field, and registered interferogram.

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Zernike term Z[4,-4] of amplitude 1.7μm generated: reconstructed wavefront, simulated interferogram, simulated far field, and registered interferogram.

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OKO Technologies supplies breadboard systems, which can be used as a basis for building of applied adaptive optics systems. Applications include scientific instrumentation, astronomy, ophthalmology, laser optics and optical alignment systems.

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Our system includes the following components:

- Deformable mirror (MMDM with up to 79 channels or PDM with up to 109 channels);
- PC interface for the deformable mirror (40-channel USB or Ethernet driver modules);
- High-voltage amplifier units for the deformable mirror;
- “FrontSurfer” WF sensor with AO loop enabled;
- All necessary cables and documentation;
- Laptop or desktop PC (optional).

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OKO Technologies also offers a complete closed-loop adaptive optical system (AOS), which includes all necessary optical and mechanical components (see Figure 21). It can be used for realtime correction of optical aberrations and generation of precision wavefronts. The speed of control depends on the camera and computer; up to 500 frames per second is achieved with inexpensive IDS uEye camera with USB3 interface and a laptop computer.

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Complete adaptive optical system with a 37-channel PDM.

References

[1] L. Baker. C tools for scientists and engineers. McGraw-Hill, 1989.

[2] C. Paterson, I. Munro, and C. Dainty. “A low cost adaptive optics system using a membrane mirror”. In: Optics Express 6 (2000), pp. 175 – 185.

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Optical Shop Testing Set-ups

原文：Flexible Optical B.V., The Netherlands（OKO Tech）

翻译：灵动智能光学，152 5716 7514 ,info@inspiroptics.com

Optical setup for testing of optical components with the Shack-Hartmann wavefront sensor should satisfy the following conditions:

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The relay optics should re-image the plane of the object to the plane of the Hartmann mask (lenslet array).

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The scheme should scale the beam is such a way that the image of the scaled object (wavefront) is densely filled with Hartmann sub-apertures.

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The optics should allow for calibration. In the general case, it consists of separate measurement of the complete setup aberration with ideal object or a source of ideal wavefront,replacing the one to be tested.

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Typical measurement scheme for testing of transparent optics with Shack- Hartmann wavefront sensor. Here CL is a collimator; L1 and L2 are lenses with focal distances f_1 and f_2, respectively; WFS is the wavefront sensor.

A typical optical optics setup for measurement of transparent optics with Shack-Hartmann wavefront sensor is shown on Figure 1. The object under test should be placed in a collimated laser beam. A telescopic system consisting of two lenses,L1 and L2,is used to re-image the aperture of the object to the microlens array(or Hartmann mask)of the wavefront sensor. The system should scale the area to be measured to the aperture size of the wavefront sensor;the scaling factor is equal tof1/f2,where f_1and f_2 are focal distances of the lenses L1 and L2,respectively.

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Although ”FrontSurfer” is able to perform “absolute” wavefront measurement without any reference pattern, it is advisable to switch to the reference mode. It allows to get rid of aberrations of the relay optics, which are always present due to imperfection of the components and misalignment in the system. Reference pattern can be measured with the object removed from the setup.

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A typical optics setup for measurement of reflective optics is shown on Figure 2. To facilitate calibration, it is convenient to build it in a similar way to the Twyman-Green interferometer, where a beam splitter divides the incoming beam in two branches. The object is placed the in first branch behind the beam splitter, and the reference (normally a high-quality flat or spherical mirror) in the second one. The beams reflected from the object and reference are them directed to the wavefront sensor. A telescopic system similar to the one from the previous scheme (Figure 1) is used to provide conjugation between the object and WFS with proper scaling. One of the branches should always be block-the first branch for measurement of the reference pattern, and the second one for testing of the object.

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Typical measurement scheme for testing of reflective optics with Shack-Hartmann wavefront sensor. See Figure 1 for notations.

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Measurement scheme for testing of a microlens with Shack-Hartmann wavefront sensor.

An example of a relay system designed to measure the aberrations of a microlens for fiber coupling applications is shown in Figure 3. In this case the microlens has a small diameter (1 mm) compared to the Hartmann mask aperture (4 mm). A telescope formed by the microlens itself and the imaging lens is used to couple the microlens aperture to the sensor mask.The Hartmann mask is co-incident with the microlens image.

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Calibration of the setup consists of two steps:
校准包括两个步骤：

(1) Obtaining of an ideal reference.
获得理想的波前参考。

(2) Calibration of the imaging microlens.
校准成像微透镜。

In the first step we ensure that the aberration of the illuminating beam is negligible over the aperture of the microlens. In the simplest case this can be achieved by illumination of the microlens with slightly divergent wavefront produced by a remote pinhole source. If there are no optical components between the pinhole and tested microlens, we can be sure the illumination wavefront is close to ideal.

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In the second step, we place another pinhole in the focus of the microlens. By doing that we eliminate all aberration of the setup, except for the aberrations introduced by the imaging lens,to

be able to measure the aberration of the relay optics only.

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In the last step, we remove the pinhole and measure the total aberration. To obtain the aberration of the microlens, the aberration of the imaging setup measured at step 2 should be extracted from the measurement result.

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In this example we have used slightly divergent wavefront instead of collimated. This does not introduce any significant problem if the radius of curvature of wavefront is much larger than the focal length of the microlens under test. If the lens has focal length of 4mm and point source is placed at 1m from the lens, spherical aberration introduced by “wrong” illumination is negligible. Nevertheless, it is always advisable to check expected aberrations using ray-tracing code.

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4. Optical conjugation （光学共轭）

The input pupil of the WF sensor should be optically conjugated with the pupil of the optical system. Conjugated is achieved by placing the sensor input in the image of the pupil. If the sensor is positioned outsider of the pupil image, any aberrations in the pupil will cause geometrical beam distortions in the sensor input, important: wrong conjugation is the cause of numerous disappointments with WF measurement.

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Usually the conjugation is achieved with an afocal telescope, formed by two positive lenses with focal lengths F1 and F2,respectively.

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The coupling telescope reduces the transversal dimension of the light beam by a factor(F1/F2), and magnifies all ray angles with the same factor, thus preserving the geometrical of aberration Δ . For the input aberration represented by a parabola $\Delta(a)=a^2/2R$ where a is the radial coordinate and R is the radius of curvature, the F number FN is defined as R/A, where A is the semi-diameter of the input pupil. Expressing the F number in terms of aberration amplitude, we obtain $FN=A/(2\Delta)$. If the aperture is reduced M times by the coupling telescope, the F number of the beam is also reduced M times.

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Scaling of the F number imposes serious practical restrictions on the range of optical powers that can be controlled by the sensor. The practically acceptable range of F numbers is always limited by the complexity and aberrations of the telescope optics While F numbers in the range between infinity and 10 are easily achievable with single lenses and achromatic doublets, the range between 10 and 3 calls for a somewhat mere complicated optics and, finally, the range between 3 and 1 can be very expensive and also very difficult to correct for intrinsic aberrations.

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It follows from Fig. 4 that the F number of the wavefront at the telescope input should be considerably smaller than (FN/M) where FN is the F number of any of two lenses forming the coupling telescope-assuming that both components of the telescope have equal F number - and M is the telescope magnification. Finally, both the wavefront curvature and the wavefront tilts should not exceed the specifications of the sensor at the input pupil of the sensor.

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Scaling of the WF in an afocal telescope with zero WF curvature(top),and compensation for a non-zero curvature (bottom).Here $\delta l=\frac{F_1^2}{R}$.

4.1. Compensation of the curvature 曲率补偿

To secure the maximum measurement range for higher-order aberrations, the input of wavefront sensor have zero spherical term. To compensate for the sphericity of incoming wavefront, a calibrated spherical or parabolic phase shift could be added to the wavefront at the input of the telescope. This can be achieved by using a pinhole as a source of diverging in the diffraction angle 1/FN.To satisfy this condition, d_p should satisfy to d_p-FN*2.44 λ.

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Also, the WF curvature in the input pupil of the sensor can be compensated by changing the distance between the lenses of the telescope, preserving the position of the input pupil in the front focus of the first lens, and keeping the point of ray convergence in the front focus of the second lens. Under these conditions, the shift of the first lens is given by $\delta l=F_1^2/R$. It is easy to see that such a compensation does not change the pupil magnification (F_1/F_2), however it can result in additional(spherical) aberration introduced into the wavefront.

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The easiest practical way to achieve such a compensation consists the wavefront sensor, together with the second lens of the telescope, on a moving stage with movement range of at least ($\pm F_1^2/R_{min}$), where R_{min} is the shortest expected radius of curvature of the input WF. In the case, to accommodate for the possible range of the input wavefront curvatures, he physical diameter of the first lens should be increased by at least ($R_{min}+F_1$)/(R_{min}) .

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4.2. Thermal stability consideration 有关热稳定性的考虑

Some camera models, particularly those with higher maximum frame rate are known to product appreciable amount of heat. Due to thermal expansion the size of the image sensor changes with the temperature increase. Thus, pixel pitch either nominal or measured during calibration might not be accurate anymore. To illustrate the Fig. 5 (top) shows typical aberration introduced by thermal expansion after changing the temperature of CMV4000 imaging chip from 30℃ to 50℃. As can be predicted main aberration type is defocus, it is measured about one order of magnitude higher than others. The measured error of radius of curvature relates linearly to the temperature of wavefront measurements, especially in absolute mode the is advised to keep the temperature of the camera as constant as possible.

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Effect of image sensor thermal expansion on wavefront measurement. Top: spurious defocus caused by increase of the CMOS sensor temperature from 30℃ to 50℃. Bottom: observed change in distance to a point source with the change of the image sensor temperature.

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
Parameter	Value
Aperture shape	approximately circular
Mirror coating	Metal or Metal + dielectric
Aperture dimensions	15mm diameter
Number of electrodes	37 (19) (see Fig. 7.6)
Control voltages V_c	0 . . . 150 to 300 V, dependent on the mirror
Initial RMS deviation from plane	less than 0.45 μm
Main initial aberration	1.5 fringes at 630nm
Maximum deflection of the mirror center	10 μm
Package dimensions	see Fig. 7.7
Weight	140 g

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Product type	Material	Pitch, μm	Radius, mm	Focal distance, mm	Maximum size
APO-Q-P2200-F209 (633)	Fused silica	2200	95,52	209	30x30 lenses (70x70mm)
APO-GB-P1500-F24.3 (633)	BK-7	1500	12,5	24,3	6x6 lenses (10x10mm)
APO-GB-P1500-F39.3 (633)	BK-7	1500	20,4	39	6x6 lenses (10x10mm)
APO-Q-P1100-F105 (633)	Fused silica	1100	48.0	105	55x55 lenses (62x62mm)
APO-Q-P1100-F209 (633)	Fused silica	1100	95,52	209	55x55 lenses (62x62mm)
APO-Q-P1061-F158 (633)	Fused silica	1061	72,4	158	55x55 lenses (62x62mm)
APO-Q-P1000-F40 (633)	Fused silica	1000	18.0	40	60x60 lenses (61x61mm)
APO-Q-P1000-F4,64 (633)	Fused silica	1000	2,12	4,64	60x60 lenses (61x61mm)
APO-Q-P500-F19,7 (633)	Fused silica	500	9	19,7	100x100 lenses (50x50mm)
APO-Q-P500-F8,31 (633)	Fused silica	500	3,8	8,31	100x100 lenses (50x50mm)
APO-Q-P300-F40 (633)	Fused silica	300	18.0	40	200x200 lenses (61x61mm)
APO-Q-P300-F118 (633)	Fused silica	300	54	118	200x200 lenses (61x61mm)
APO-Q-P300-F2 (633)	Fused silica	300	0,94	2,0	200x200 lenses (61x61mm)
APO-Q-P384-F5.75 (633)	Fused silica	384	2 627	5,75	5x5 lenses (2,5x2,5mm)
APO-Q-P222-F1,86 (633)	Fused silica	222	0,85	1,86	200x200 lenses (45x45mm)
APO-Q-P222-F0,93 (633)	Fused silica	222	0,43	0,93	200x200 lenses (45x45mm)
APO-Q-P192-F5.75 (633)	Fused silica	192	2 627	5,75	9x9 lenses (2x2mm)
APO-Q-P192-F3.17 (633)	Fused silica	192	1,45	3,17	200x200 lenses (45x45mm)
APO-Q-P200-F7.0 (633)	Fused silica	200	3,2	7.0	300x300 lenses (61x61mm)
APO-Q-P200-F40 (633)	Fused silica	200	18,2	40	300x300 lenses (61x61mm)
APO-Q-P150-F10 (633)	Fused silica	150	4,57	10	400x400 lenses (61x61mm)
APO-Q-P150-F3.5 (633)	Fused silica	150	1,6	3,5	400x400 lenses (61x61mm)
APO-Q-P148-F1,24 (633)	Fused silica	148	0,6	1,24	400x400 lenses (61x61mm)
APO-Q-P100-F0,217 (532)	Fused silica	100	0,1	0,217	500x500 lenses (51x51mm)
APO-GT-P100-F0,143 (633)	StiH-53	100	0,12	0,143	500x500 lenses (51x51mm)
APO-P(GB)-P250-F0,57 (633)	BK-7/NOA-61	250	0,32	0,57	25x25mm
APO-P(GB)-P200-F1,25 (633)	BK-7/NOA-61	200	0,7	1,25	40x40mm
APO-P(GB)-P200-F6,4 (633)	BK-7/NOA-61	200	3,32	6,4	25x25mm
APO-P(GB)-P300-F1,7 (633)	BK-7/NOA-61	300	0,94	1,7	25x25mm
APO-P(GB)-P300-F107 (633)	BK-7/NOA-61	300	60	107	25x25mm

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- Antonello, Jacopo, et al. "Optimization-based wavefront sensorless adaptive optics for multiphoton microscopy." JOSA A 31.6 (2014): 1337-1347.
- Van Werkhoven, T. I. M., et al. "Snapshot coherence-gated direct wavefront sensing for multi-photon microscopy." Optics express 22.8 (2014): 9715-9733.
- Popovic, Zoran, et al. "Dual conjugate adaptive optics prototype for wide field high resolution retinal imaging." Adaptive Optics Progress. InTechOpen, 2012.
- Niu, Saisai, et al. "High-resolution retinal imaging with micro adaptive optics system." Applied optics 50.22 (2011): 4365-4375.
- Bortoletto, Favio, et al. "Multiphoton fluorescence microscopy with GRIN objective aberration correction by low order adaptive optics." PloS one 6.7 (2011): e22321.
- Thaung, Jörgen, et al. "Dual-conjugate adaptive optics for wide-field high-resolution retinal imaging." Optics express 17.6 (2009): 4454-4467.
- Wright, A. J., et al. "Adaptive optics for enhanced signal in CARS microscopy." Optics express 15.26 (2007): 18209-18219.
- Merino, David, et al. "Adaptive optics enhanced simultaneous en-face optical coherence tomography and scanning laser ophthalmoscopy." Optics express 14.8 (2006): 3345-3353.
- Wright, Amanda J., et al. "Dynamic closed-loop system for focus tracking using a spatial light modulator and a deformable membrane mirror." Optics express 14.1 (2006): 222-228.
- Rueckel, Markus, Julia A. Mack-Bucher, and Winfried Denk. "Adaptive wavefront correction in two-photon microscopy using coherence-gated wavefront sensing." Proceedings of the National Academy of Sciences 103.46 (2006): 17137-17142.
- Albert, O., et al. "Smart microscope: an adaptive optics learning system for aberration correction in multiphoton confocal microscopy." Optics letters 25.1 (2000): 52-54.

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- Cheung, D. C. L., T. H. Barnes, and T. G. Haskell. "Feedback interferometry with membrane mirror for adaptive optics." Optics communications 218.1-3 (2003): 33-41.
- Czarske, Jürgen W., et al. "Smart laser interferometer with electrically tunable lenses for flow velocity measurements through disturbing interfaces." Photonics. Vol. 2. No. 1. Multidisciplinary Digital Publishing Institute, 2015.
- Büttner, Lars, Christoph Leithold, and Jürgen Czarske. "Advancement of an Interferometric Flow Velocity Measurement Technique by Adaptive Optics." International Journal of Optomechatronics 8.1 (2014): 1-13.
- Büttner, Lars, Christoph Leithold, and Jürgen Czarske. "Interferometric velocity measurements through a fluctuating gas-liquid interface employing adaptive optics." Optics express 21.25 (2013): 30653-30663.

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- Vdovin, Gleb, and Vadim Kiyko. "Intracavity control of a 200-W continuous-wave Nd: YAG laser by a micromachined deformable mirror." Optics Letters 26.11 (2001): 798-800.
- Zeek, Erik, et al. "Adaptive pulse compression for transform-limited 15-fs high-energy pulse generation." Optics letters 25.8 (2000): 587-589.
- Hofmann, Oskar, et al. "Model-based analysis of highly dynamic laser beam shaping using deformable mirrors." Procedia CIRP 74 (2018): 602-606.
- Maxson, Jared, et al. "Adaptive electron beam shaping using a photoemission gun and spatial light modulator." Physical Review Special Topics-Accelerators and Beams 18.2 (2015): 023401.
- Garduño-Mejía, Jesus, Alan H. Greenaway, and Derryck T. Reid. "Programmable spectral phase control of femtosecond pulses by use of adaptive optics and real-time pulse measurement." JOSA B 21.4 (2004): 833-843.
- Garduño-Mejía, J., A. H. Greenaway, and D. T. Reid. "Designer femtosecond pulses using adaptive optics." Optics Express 11.17 (2003): 2030-2040.
- Armstrong, Michael R., et al. "Versatile 7-fs optical parametric pulse generation and compression by use of adaptive optics." Optics Letters 26.15 (2001): 1152-1154.
- Bartels, Randy, et al. "Shaped-pulse optimization of coherent emission of high-harmonic soft X-rays." Nature 406.6792 (2000): 164-166.
- Zeek, Erik, et al. "Pulse compression by use of deformable mirrors." Optics letters 24.7 (1999): 493-495.
- Druon, Frédéric, et al. "Wave-front correction of femtosecond terawatt lasers by deformable mirrors." Optics letters 23.13 (1998): 1043-1045.

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- Loktev, Mikhail, et al. "Adaptive optics combined with computer post-processing for horizontal turbulent imaging." Real-Time Image and Video Processing 2012. Vol. 8437. International Society for Optics and Photonics, 2012.
- Loktev, Mikhail, et al. "Speckle imaging through turbulent atmosphere based on adaptable pupil segmentation." Optics letters 36.14 (2011): 2656-2658.
- Vdovin, Gleb, et al. "Imaging through turbulence with temporally and spatially multiplexed systems." Optics in Atmospheric Propagation and Adaptive Systems XIV. Vol. 8178. International Society for Optics and Photonics, 2011.
- Loktev, M., et al. "Experiments on speckle imaging using projection methods." Unconventional Imaging, Wavefront Sensing, and Adaptive Coded Aperture Imaging and Non-Imaging Sensor Systems. Vol. 8165. International Society for Optics and Photonics, 2011.
- Keller, Christoph U., Claude Plymate, and S. Mark Ammons. "Low-cost solar adaptive optics in the infrared." Innovative Telescopes and Instrumentation for Solar Astrophysics. Vol. 4853. International Society for Optics and Photonics, 2003.

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- Zhou G. Unveiling coronal fine structures with adaptive optics: Solar physics[J]. Nature Astronomy, 2025, 9(8): 1105-1106.
- Fan R, Wang Z, Li P, et al. Metasurfaces in adaptive optics: a new opportunity in optical wavefront sensing[J]. Laser & Photonics Reviews, 2025, 19(24): e01566.
- Zhou Y, Li L, Zhang J, et al. Meta‐Device for Field‐of‐View Tunability via Adaptive Optical Spatial Differentiation[J]. Advanced Science, 2025, 12(9): 2412794.
- Hou W, He Y, Shen Y, et al. Deep Learning‐Based Computational Adaptive Optics for Photoacoustic Microscopy[J]. Laser & Photonics Reviews, 2026, 20(4): e01943.
- Tobón-Maya H, Willstatter L, Zapata-Valencia S I, et al. Hadamard based single-pixel microscopy using sensor-less adaptive optics supported by multi-actuator adaptive lens[J]. Nature Communications, 2025, 16(1): 10933.
- Cabillic M, Forriere H, Bettarel L, et al. In-depth single molecule localization microscopy using adaptive optics and single objective light-sheet microscopy[J]. Nature Communications, 2025, 16(1): 8362.
- Rao C, Zhong L, Guo Y, et al. Astronomical adaptive optics: a review[J]. PhotoniX, 2024, 5(1): 16.
- Liu Y, Crowell J A, Kurokawa K, et al. Ultrafast adaptive optics for imaging the living human eye[J]. Nature Communications, 2024, 15(1): 10409.
- Streich, Lina, et al. "High-resolution structural and functional deep brain imaging using adaptive optics three-photon microscopy." Nature methods 18 (2021): 1253‐1258.
- Hao, Xiang, et al. "Three-dimensional adaptive optical nanoscopy for thick specimen imaging at sub-50-nm resolution." Nature methods 18.6 (2021): 688-693.
- Lin, Ruizhe, et al. "Subcellular three-dimensional imaging deep through multicellular thick samples by structured illumination microscopy and adaptive optics." Nature communications 12.1 (2021): 1-14.
- Rodríguez, Cristina, et al. "An adaptive optics module for deep tissue multiphoton imaging in vivo." Nature Methods 18.10 (2021): 1259-1264.
- Lu, Rongwen, et al. "In-vivo sub-diffraction adaptive optics imaging of photoreceptors in the human eye with annular pupil illumination and sub-Airy detection." Optica 8.3 (2021): 333-343.
- Wu, Jiamin, et al. "Iterative tomography with digital adaptive optics permits hour-long intravital observation of 3D subcellular dynamics at millisecond scale." Cell (2021).
- Zhou Zhou,et al.” Adaptive optical microscopy via virtual imaging assisted wavefront sensing for high resolution tissue imaging”, PhotoniX (2022) 3:13 1-20